

Ez SHOT HD





	Good Doctors Wireless Solution	Caries Detector	Wireless Intraoral Camera	Wired Intraoral Camera	Apex Locator	Curing Light
	IC-AD	IC-WHCD100 IC-DCDU75	IC-WHT60 IC-DCAM70 IC-EZAM80	IC-WHMA90 IC-EZAU81 IC-WHW61 IC-DCCU74 IC-DCAU71 IC-EZBU82	AL-DFA20 AL-DF10	CL-DP40 CL-DC20 CL-DP41 CL-DL10
		А		O USE intraoral camera	a	
	eatures of Ez SHO	r hd		_		
Wire	reless type only High speed digital wireless High speed digital wireless provide more clear live image without interruption.					- 1
	Wireless receiver Compact wireless receiver without any power supply.					
******	Macro to Full arch Camera is able to focus from	n macro mode (close) to full arch.			
	Auto focus Can get clear image with auto focus function in oral cavity.					
FHD	FHD image support The device shows FHD (1080p) in live image supported by a 5M sensor.					
\$**** 1 1 \$_{k_{2}}\$*	Shadowless Lighting Rounded 12pcs LEDs make shadowless lighting.					
×	Adjustable LED brightness Image quality loss by light reflection, can adjust LED brightness.					
	Black & White function Can take and send analog X-ray film image to PC.					
	Plug & Play					T

Plug & Play No need additional driver to install the device.

Easy to capture

Easy to image capture by camera only without extra device like foot pedal.

Software compatibility

Can be use camera with most of dental imaging software. * It would need twain driver for some of dental imaging software.

2 Specification

ľ

Model name		IC-EZAM80(Wireless model)
Camera	Operation voltage	3.7V
	Size	Approx. 255 x 27 x 29 mm
	Weight	Approx. 89g
	Resolution	FHD (1920×1080)
	Sensor	5M Sensor CMOS / Liquid
Charger	Operation voltage	5V
	Size	Approx. 110 x 105 x 68 mm
	Weight	Approx. 100g
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

Model name		IC-EZAU81(Wired model)	
Camera	Operation voltage	5V	
	Size	Approx. 202 x 27 x 29 mm	
	Weight	Approx. 55g	
	Resolution	FHD (1920x1080)	
	Sensor	5M Sensor CMOS / Liquid	
USB cable	Length	Approx. 3M	
Disposable cover	Size	Approx. 32 x 100 mm	
	Quantity	Approx. 100 pcs by weight x 2ea	
	Raw material	Polyethylene	

3 Minimum System Requirement

PC system:	3.3 Ghz CPU, 4GB RAM	
OS:	Windows 7	
USB port:	USB 2.0	



Upper images are take by real product.

Wired



GoodArs

www.gooddrs.com E-mail. info@gooddrs.com



Good Doctors Co., Ltd. #208, B-dong, 283 Bupyeong-daero, Bupyeong-gu (Woolim Lions Valley, Cheongchoen-dong), Incheon, 21315, Republic of Korea Tel: +82 32 424 6325 Fax: +82 32 424 6326

Wireless